

Ctrb1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P17538

Item Number: CM0016097

Introduction

CTRB1 Rabbit Polyclonal Antibody (pAb) validated for WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, rat Chymotrypsinogen B (CTRB1). Suitable for research on serine protease secretion and digestive enzyme biology.

[Learn More](#)

Parameter	Value
Product Name	CTRB1 Rabbit pAb
Remarks/Alias	CTRB; CTRB1
Species	Human
GeneID (human)	1504
GeneID	1504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 79-263 of human CTRB1 (NP_001897.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17538
Protein Weight	28kDa
Shipping	Ice bag