

Cxcr3-B Rabbit Pab - P49682-2

Item Number: CM0024302

Introduction

CXCR3-B Rabbit Polyclonal Antibody for WB, IF-P, ELISA. Recognizes human and mouse CXCR3-B isoform 2, a chemokine receptor inhibiting angiogenesis via CXCL4/CXCL10. Useful for oncology and vascular studies.

[Learn More](#)

Parameter	Value
Product Name	CXCR3-B Rabbit pAb
Remarks/Alias	GPR9; MigR; CD182; CD183; Mig-R; CKR-L2; CMKAR3; IP10-R; CXCR3
Species	Human
Gene ID (Human)	2833
Gene ID	2833
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 1-100 human CXCR3-B(NP_001136269.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49682-2
Protein Weight	45kDa
Shipping	Ice bag