


Cxcr7 Rabbit Mab For Wb, Elisa - P25106

Item Number: CM0005923



Introduction

This CXCR7 rabbit monoclonal antibody (WB, ELISA) targets human atypical chemokine receptor 3, suitable for cancer and neuroscience research. Detects ~41 kDa band; validated in human samples.

[Learn More](#)

Parameter	Specification
Product Name	CXCR7 Rabbit mAb
Remarks/Alias	RDC1; CXCR7; RDC-1; CMKOR1; CXC-R7; CXCR-7; GPR159
Species	Human
GeneID (Human)	57007
GeneID	57007
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 263-362 of human CXCR7 (NP_064707.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P25106
Protein Weight	41kDa
Shipping	Ice bag