

Cy3 Rabbit Monoclonal Antibody - Cy3

Item Number: CM0005647

Introduction

This rabbit monoclonal antibody targets the Cy3 cyanine dye (CAS 146397-20-8), offering species-independent reactivity for DB and ELISA assays. Suitable for direct labeling or anti-dye signal generation strategies.

[Learn More](#)

Parameter	Value
Product Name	Cy3 Rabbit mAb
Remarks / Alias	Cy3
Gene ID (Human)	CAS:146397-20-8
Gene ID	CAS:146397-20-8
Immunogen	Chemical compounds corresponding to Cy3
Source	Rabbit
Category	Monoclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag