

Cy3/Cy5 Rabbit Mab

Item Number: CM0005240

Introduction

Cy3/Cy5 Rabbit monoclonal antibody for DB and ELISA applications. Exhibits species-independent cross-reactivity, suitable for universal detection of Cy3/Cy5 conjugates. Sourced from rabbit, ideal for immunoassay development and diagnostic raw material supply.

[Learn More](#)

Parameter	Specification
Product Name	Cy3/Cy5 Rabbit mAb
Gene ID (CAS)	CAS:146368-15-2
Immunogen	Chemical compounds corresponding to Cy5
Source	Rabbit
Category	Monoclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag