


Cyclin B1 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - P14635

Item Number: CM0006441



Introduction

High-quality Cyclin B1 rabbit monoclonal antibody validated for WB, IP, IHC-P, ELISA. Reacts with human Cyclin B1 (P14635), essential for G2/M cell cycle transition. Ideal for cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	Cyclin B1 Rabbit mAb
Remarks/Alias	CCNB; Cyclin B1
Species	Human
GeneID (Human)	891
GeneID	891
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 334-433 of human Cyclin B1 (P14635).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P14635
Protein Weight	48 kDa
Shipping	Ice bag