


Cyclin B1 Rabbit Pab For Wb, Ihc-P, If/Icc, Ip, Elisa - P14635

Item Number: CM0020889



Introduction

Cyclin B1 Rabbit pAb detects endogenous Cyclin B1 protein in human, mouse, rat samples. Suitable for WB, IHC-P, IF/ICC, IP, and ELISA. Target is a key regulator of G2/M transition and mitosis.

[Learn More](#)

Parameter	Value
Product Name	Cyclin B1 Rabbit pAb
Remarks/Alias	CCNB; Cyclin B1
Species	Human
GeneID (Human)	891
GeneID	891
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 334-433 of human Cyclin B1 (NP_114172.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14635
Protein Weight	48kDa
Shipping	Ice bag