


## Cyclin D3 Rabbit Mab - P30281

Item Number: CM0005277



### Introduction

Recombinant rabbit monoclonal antibody targeting human Cyclin D3 (CCND3), validated for WB and ELISA, recognizing human, mouse, rat. Ideal for cell cycle research.

[Learn More](#)

| Parameter        | Value                       |
|------------------|-----------------------------|
| Product Name     | Cyclin D3 Rabbit mAb        |
| Remarks/Alias    | CCND3; cyclin D3; Cyclin D3 |
| Species          | Human                       |
| GeneID (Human)   | 896                         |
| GeneID           | 896                         |
| Immunogen        | Recombinant protein         |
| Source           | Rabbit                      |
| Category         | Monoclonal Antibodies       |
| Application      | WB, ELISA                   |
| Cross Reactivity | Human, Mouse, Rat           |
| SWISS            | P30281                      |
| Protein Weight   | 33kDa                       |
| Shipping         | Ice bag                     |