

Cyclin E1 Rabbit Polymab® For Wb And Elisa - P24864

Item Number: CM0006879



Introduction

Rabbit monoclonal antibody targeting human Cyclin E1 (CCNE1), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for cell cycle and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	Cyclin E1 Rabbit PolymAb®
Remarks/Alias	CCNE; pCCNE1
Species	Human
Gene ID (Human)	898
Gene ID	898
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human Cyclin E1 (NP_001229.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24864
Protein Weight	47 kDa
Shipping	Ice bag