

Cyp11A1 Rabbit Polyclonal Antibody (Wb, If, Ihc, Elisa) - P05108

Item Number: CM0022531



Introduction

Rabbit polyclonal antibody against CYP11A1 (P450scc), validated for Western blot, immunofluorescence, IHC and ELISA. Reacts with human, mouse, and rat CYP11A1, a mitochondrial cholesterol side-chain cleavage enzyme essential for steroidogenesis.

[Learn More](#)

Parameter	Value
Product Name	CYP11A1 Rabbit pAb
Remarks/Alias	CYP11A; CYPXIA1; P450SCC; CYP11A1
Species	Human
GeneID	1583
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05108
Protein Weight	60kDa
Shipping	Ice bag