

Cyp24A1 Rabbit Monoclonal Antibody For IHC-P, IF/ICC, Elisa - Q07973

Item Number: CM0003087

Introduction

CYP24A1 Rabbit Monoclonal Antibody validated for IHC-P, IF/ICC, and ELISA. Targets human CYP24A1 (Q07973), a mitochondrial enzyme critical for vitamin D catabolism and calcium homeostasis. Research use only.

[Learn More](#)

Parameter	Value
Product Name	CYP24A1 Rabbit mAb
Remarks/Alias	CP24, HCAI, CYP24, HCINF1, P450-CC24, CYP24A1
Species	Human
Gene ID (Human)	1591
Gene ID	1591
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-410 of human CYP24A1 (NP_000773.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q07973
Protein Weight	59kDa
Shipping	Ice bag