

Cyp24A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q07973

Item Number: CM0019745

Introduction

Rabbit polyclonal CYP24A1 antibody validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat. Suitable for vitamin D catabolism, calcium homeostasis, and mitochondrial 24-hydroxylase function studies.

[Learn More](#)

Parameter	Value
Product Name	CYP24A1 Rabbit pAb
Remarks/Alias	CP24; HCAI; CYP24; HCINF1; P450-CC24; CYP24A1
Species	Human
Gene ID (Human)	1591
Gene ID	1591
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 36-448 of human CYP24A1 (NP_001122387.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07973
Protein Weight	59kDa
Shipping	Ice bag