

Cyp7A1 Mouse Monoclonal Antibody For Wb, Elisa - P22680

Item Number: CM0006540



Introduction

Mouse monoclonal antibody against human CYP7A1, validated for WB and ELISA. Cross-reacts with mouse and rat. Recombinant immunogen (aa 86-108). For cholesterol metabolism and bile acid synthesis research.

[Learn More](#)

Parameter	Value
Product Name	CYP7A1 Mouse mAb
Remarks/Alias	CP7A; CYP7; CYPVII; CYP7A1
Species	Human
GeneID	1581 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 86-108 of human CYP7A1 (P22680)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P22680
Protein Weight	58 kDa
Shipping	Ice bag