


Cystatin C Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P01034

Item Number: CM0014675



Introduction

Cystatin C rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects endogenous Cystatin C in human, mouse, rat samples. Ideal for renal function, amyloidosis, and protease inhibitor research.

[Learn More](#)

Specification	Details
Product Name	Cystatin C Rabbit pAb
Remarks/Alias	ARMD11; HEL-S-2; Cystatin C
Species	Human
Gene ID	1471 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human Cystatin C (NP_000090.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01034
Protein Weight	13kDa
Shipping	Ice bag