

Cytochrome C Rabbit Mab - P99999

Item Number: CM0011288



Introduction

Anti-Cytochrome c rabbit monoclonal antibody for WB, IHC-P, ELISA. Detects Cytochrome c (CYC, HCS, THC4) across Human, Mouse, Rat. Molecular weight: 12kDa. Ideal for apoptosis and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	Cytochrome c Rabbit mAb
Remarks/Alias	CYC; HCS; THC4
Species	Human
GeneID (human)	54205
GeneID	54205
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human Cytochrome C (NP_061820.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P99999
Protein Weight	12kDa
Shipping	Ice bag