


Cytokeratin 14 (Krt14) Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P02533

Item Number: CM0011401



Introduction

Rabbit monoclonal antibody targeting human Cytokeratin 14 (KRT14). Validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for epithelial biology and basal keratinocyte research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 14 (KRT14) Rabbit mAb
Remarks/Alias	K14; NFJ; CK14; EBS1; EBS3; EBS4; EBS1A; EBS1B; EBS1C; EBS1D; Cytokeratin 14 (KRT14)
Species	Human
GeneID (Human)	3861
GeneID	3861
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 373-472 of human Cytokeratin 14 (KRT14)(NP_000517.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02533
Protein Weight	52kDa
Shipping	Ice bag