

Cytokeratin 17 (Krt17) Rabbit Monoclonal Antibody For Wb, Elisa - Q04695

Item Number: CM0005997

Introduction

Validated for WB and ELISA, this rabbit monoclonal antibody targets human Cytokeratin 17 (KRT17), cross-reacting with mouse. KRT17 is a type I keratin essential for skin appendage formation and a marker in psoriasis and epithelial cancers.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 17 (KRT17) Rabbit mAb
Remarks/Alias	PC; K17; PC2; 39.1; CK-17; PCHC1; Cytokeratin 17 (KRT17)
Species	Human
GeneID	3872
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 332-432 of human Cytokeratin 17 (KRT17). (NP_000413.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q04695
Protein Weight	48 kDa
Shipping	Ice bag