

Cytokeratin 17 (Krt17) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q04695

Item Number: CM0014108

Introduction

High-quality rabbit polyclonal antibody against human Cytokeratin 17 (KRT17) for use in Western blot (WB), immunofluorescence/immunocytochemistry (IF/ICC), and ELISA. Cross-reacts with mouse and rat. Ideal for epithelial biology and dermatology research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 17 (KRT17) Rabbit pAb
Remarks/Alias	PC; K17; PC2; 39.1; CK-17; PCHC1; Cytokeratin 17 (KRT17)
Species	Human
Gene ID (Human)	3872
Gene ID	3872
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 133-432 of human Cytokeratin 17 (KRT17) (NP_000413.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04695
Protein Weight	48kDa
Shipping	Ice bag