

Cytokeratin 18 (Krt18) Rabbit Monoclonal Antibody For Wb - P05783

Item Number: CM0003612

Introduction

Cytokeratin 18 (KRT18) rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat KRT18. Ideal for epithelial cell biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 18 (KRT18) Rabbit mAb
Remarks/Alias	K18; CK-18; CYK18; Cytokeratin 18 (KRT18)
Species	Human
GenelD (Human)	3875
GenelD	3875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-430 of human Cytokeratin 18 (KRT18) (P05783)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05783
Protein Weight	48kDa
Shipping	Ice bag