

Cytokeratin 20 (Ck20) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P35900

Item Number: CM0002098

Introduction

Cytokeratin 20 (CK20) Rabbit mAb validated for WB, IHC-P, ELISA. Recognizes human and rat CK20 (~48 kDa), expressed in intestinal epithelia and Merkel cells. Ideal for gastrointestinal tract research and IHC marker studies.

[Learn More](#)

Parameter	Specification
Product Name	Cytokeratin 20 (CK20) Rabbit mAb
Remarks/Alias	K20; CD20; KRT20; CK-20; KRT21; Cytokeratin 20 (CK20)
Species	Human
Gene ID (Human)	54474
Gene ID	54474
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P35900
Protein Weight	48 kDa
Shipping	Ice bag