

Dab2Ip Rabbit Pab - Q5Vwq8

Item Number: CM0019581

Introduction

High-quality rabbit polyclonal antibody against DAB2IP, a scaffold protein involved in signaling pathways, apoptosis, and tumor suppression. Validated for WB, IF/ICC, ELISA in human and rat samples.

[Learn More](#)

Parameter	Value
Product Name	DAB2IP Rabbit pAb
Remarks/Alias	AIP1; AIP-1; AF9Q34; DIP1/2; DAB2IP
Species	Human
GeneID (Human)	153090
GeneID	153090
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 205-320 of human DAB2IP (NP_619723.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q5VWQ8
Protein Weight	132kDa
Shipping	Ice bag