

Dap Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51397

Item Number: CM0030232

Introduction

Explore CamelBio's DAP rabbit polyclonal antibody (P51397) for Western blot, IHC-P, IF/ICC, and ELISA. Validated in human, mouse, rat. Ideal for apoptosis, autophagy, and translation studies.

[Learn More](#)

Parameter	Value
Product Name	DAP Rabbit pAb
Remarks/Alias	DAP1; DAP
Species	Human
GeneID (Human)	1611
GeneID	1611
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-102 of human DAP (NP_004385.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51397
Protein Weight	11kDa
Shipping	Ice bag