

DAZL Rabbit Pab - Q92904

Item Number: CM0013307

Introduction

DAZL Rabbit Polyclonal Antibody validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. RNA-binding protein essential for spermatogenesis; testis-specific expression. Ideal for reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	DAZL Rabbit pAb
Remarks/Alias	DAZH; DAZL1; DAZLA; SPGYLA; DAZL
Species	Human
GenelD (Human)	1618
GenelD	1618
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92904
Protein Weight	33kDa
Shipping	Ice bag