

Dcaf11 Rabbit Pab - Q8Teb1

Item Number: CM0021607

Introduction

Rabbit polyclonal antibody against DCAF11, a substrate receptor for CUL4-DDB1 E3 ubiquitin ligase complex. Validated for WB and ELISA in human, mouse, rat. Useful for ubiquitination pathway research.

[Learn More](#)

Parameter	Value
Product Name	DCAF11 Rabbit pAb
Remarks/Alias	GL014; WDR23; PRO2389; DCAF11
Species	Human
GeneID (Human)	80344
GeneID	80344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human DCAF11 (NP_079506.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TEB1
Protein Weight	62kDa
Shipping	Ice bag