

Ddx1 Rabbit Mab For Wb, Ihc-P, Ip, Elisa - Q92499

Item Number: CM0010454



Introduction

Anti-DDX1 rabbit monoclonal antibody validated for WB, IHC-P, IP, ELISA with cross-reactivity to human, mouse, rat. Ideal for studying RNA helicase roles in transcription, antiviral response, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	DDX1 Rabbit mAb
Remarks/Alias	DBP-RB; UKVH5d; DDX1
Species	Human
GeneID (Human)	1653
GeneID	1653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 441-740 of human DDX1 (NP_004930.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92499
Protein Weight	82 kDa
Shipping	Ice bag