


Ddx3 Rabbit Pab For Wb, If/Icc, Elisa - 000571

Item Number: CM0030384



Introduction

High-quality rabbit polyclonal antibody targeting human DDX3 (000571) for Western blot, IF/ICC, and ELISA. Involved in RNA metabolism, translation regulation, and innate immune signaling.

[Learn More](#)

Parameter	Value
Product Name	DDX3 Rabbit pAb
Remarks/Alias	DBX; DDX3; HLP2; DDX14; CAP-Rf; MRX102; MRXSSB; DDX3X
Species	Human
Gene ID (Human)	1654
Gene ID	1654
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human DDX3 (NP_001347.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	000571
Protein Weight	73kDa
Shipping	Ice bag