


Ddx39A Rabbit Polyclonal Antibody For Wb And Elisa - 000148

Item Number: CM0012823



Introduction

DDX39A Rabbit Polyclonal Antibody validated for Western blotting and ELISA; cross-reacts with mouse and rat. Target ATP-dependent RNA helicase involved in mRNA export and splicing. Ideal for nucleocytoplasmic transport and RNA metabolism research.

[Learn More](#)

Parameter	Details
Product Name	DDX39A Rabbit pAb
Remarks/Alias	BAT1; DDXL; BAT1L; DDX39; URH49; DDX39A
Species	Human
GeneID (Human)	10212
GeneID	10212
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 178-427 of human DDX39A (NP_005795.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O00148
Protein Weight	49kDa
Shipping	Ice bag