


Ddx50 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bq39

Item Number: CM0025593



Introduction

Rabbit polyclonal antibody targeting human DDX50 (GU2/GUB), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for RNA metabolism and antiviral restriction research.

[Learn More](#)

Parameter	Value
Product Name	DDX50 Rabbit pAb
Remarks/Alias	GU2; GUB; mcdrh; RH-II/GuB; DDX50
Species	Human
GeneID (Human)	79009
GeneID	79009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human DDX50 (NP_076950.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQ39
Protein Weight	83kDa
Shipping	Ice bag