

Ddx51 Rabbit Pab - Q8N8A6

Item Number: CM0016080

Introduction

DDX51 Rabbit Polyclonal Antibody validated for WB and ELISA with cross-reactivity across Human, Mouse, and Rat. Target: ATP-dependent RNA helicase DDX51, a key factor in 60S ribosomal subunit biogenesis.

[Learn More](#)

Parameter	Value
Product Name	DDX51 Rabbit pAb
Remarks/Alias	DDX51
Species	Human
GeneID (human)	317781
GeneID	317781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 367-666 of human DDX51 (NP_778236.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N8A6
Protein Weight	72kDa
Shipping	Ice bag