

Ddx55 Rabbit Pab For Wb, Elisa - Q8Nhq9

Item Number: CM0018961

Introduction

DDX55 rabbit polyclonal antibody for Western blot and ELISA. Detects human and mouse DDX55, a DEAD-box RNA helicase involved in rRNA processing. Ideal for ribosome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	DDX55 Rabbit pAb
Remarks/Alias	DDX55
Species	Human
GeneID (human)	57696
GeneID	57696
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-360 of human DDX55 (NP_065987.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NHQ9
Protein Weight	69kDa
Shipping	Ice bag