


Decr1 Rabbit Pab - Q16698

Item Number: CM0012845



Introduction

DEC1 Rabbit pAb is a polyclonal antibody for WB and ELISA, targeting human, mouse, rat DEC1 (mitochondrial 2,4-dienoyl-CoA reductase). Ideal for fatty acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	DEC1 Rabbit pAb
Remarks/Alias	DEC1; NADPH; SDR18C1; DEC1
Species	Human
GeneID (Human)	1666
GeneID	1666
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 35-335 of human DEC1 (NP_001350.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16698
Protein Weight	36kDa
Shipping	Ice bag