


Decr2 Rabbit Polyclonal Antibody (Pab) For Wb, Elisa - Q9Nui1

Item Number: CM0019133



Introduction

DEC2 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Detects a 31 kDa band. Targets peroxisomal 2,4-dienoyl-CoA reductase, an auxiliary enzyme in fatty acid β -oxidation.

[Learn More](#)

Parameter	Value
Product Name	DEC2 Rabbit pAb
Remarks/Alias	PDCR; SDR17C1; DEC2
Species	Human
Gene ID (Human)	26063
Gene ID	26063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-240 of human DEC2 (NP_065715.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NUI1
Protein Weight	31kDa
Shipping	Ice bag