

Defb4A Rabbit Pab For Ihc-P, Elisa - O15263

Item Number: CM0018872



Introduction

DEFB4A Rabbit pAb is a rabbit polyclonal antibody targeting human beta-defensin 2, validated for IHC-P and ELISA with cross-reactivity to mouse and rat. Used in antimicrobial peptide research and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	DEFB4A Rabbit pAb
Remarks/Alias	BD-2; SAP1; DEFB2; DEFB4; HBD-2; DEFB-2; DEFB102; DEFB4A
Species	Human
GeneID (Human)	1673
GeneID	1673
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-64 of human DEFB4A (NP_004933.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15263
Protein Weight	7kDa
Shipping	Ice bag