


Dek Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - P35659

Item Number: CM0002369



Introduction

DEK Rabbit monoclonal antibody for WB, IHC-P, IF/ICC, ELISA in human, rat. Targets nuclear chromatin regulator DEK (43kDa). Ideal for cancer and chromatin research. - P35659

[Learn More](#)

Parameter	Value
Product Name	DEK Rabbit mAb
Remarks/Alias	D6S231E; DEK
Species	Human
GeneID (Human)	7913
GeneID	7913
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-200 of human DEK (P35659).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P35659
Protein Weight	43kDa
Shipping	Ice bag