


Depdc6 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tb45

Item Number: CM0025270



Introduction

High-quality DEPDC6 rabbit polyclonal antibody validated for Western blot (WB) and ELISA applications. Detects human and mouse DEPTOR, a key negative mTOR regulator essential for autophagy studies. Ideal for IVD research and development.

[Learn More](#)

Parameter	Value
Product Name	DEPDC6 Rabbit pAb
Remarks/Alias	DEP.6; DEPDC6
Species	Human
GeneID (Human)	64798
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-217 of human DEPDC6 (NP_073620.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TB45
Protein Weight	46kDa
Shipping	Ice bag