

Desmin Rabbit Polyclonal Antibody For Wb, Elisa - P17661

Item Number: CM0013865

Introduction

Desmin rabbit polyclonal antibody validated for WB and ELISA, cross-reacts with mouse. Targets human Desmin, a type III intermediate filament essential for muscle sarcomere integrity. Suitable for cardiac and skeletal muscle research.

[Learn More](#)

Parameter	Value
Product Name	Desmin Rabbit pAb
Remarks/Alias	CSM1; CSM2; CDCD3; LGMD1D; LGMD1E; LGMD2R; Desmin
Species	Human
GeneID (Human)	1674
GeneID	1674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Desmin (NP_001918.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P17661
Protein Weight	54kDa
Shipping	Ice bag