

Desmocollin 3 (Dsc3) Rabbit Polyclonal Antibody For Wb, Elisa - Q14574

Item Number: CM0022180

Introduction

Rabbit polyclonal antibody targeting human Desmocollin 3 (DSC3), designed for Western blot and ELISA applications. Cross-reacts with rat. Suitable for studying desmosome adhesion and epidermal barrier integrity.

[Learn More](#)

Parameter	Value
Product Name	Desmocollin 3 Rabbit pAb
Remarks/Alias	DSC; DSC1; DSC2; DSC4; CDHF3; HT-CP; Desmocollin 3
Species	Human
Gene ID (Human)	1825
Gene ID	1825
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 576-690 of human Desmocollin 3 (NP_001932.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q14574
Protein Weight	100kDa
Shipping	Ice bag