


Desmoplakin Rabbit Mab For Wb, If-P, Ihc-P, Elisa - P15924

Item Number: CM0007972



Introduction

Desmoplakin Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat. Targets the ~332 kDa desmosomal protein, ideal for cell adhesion, cardiomyopathy, and dermatology research.

[Learn More](#)

Parameter	Value
Product Name	Desmoplakin Rabbit mAb
Remarks/Alias	DP; DCWHKTA
Species	Human
GeneID (Human)	1832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2609-2822 of human Desmoplakin (NP_004406.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15924
Protein Weight	332kDa
Shipping	Ice bag