

Destrin Rabbit Monoclonal Antibody For Wb And Elisa - P60981

Item Number: CM0003481

Introduction

Destrin rabbit monoclonal antibody validated for WB and ELISA applications, cross-reactive with human, mouse, and rat. Suitable for investigating actin depolymerization and cytoskeletal regulation. UniProt P60981.

[Learn More](#)

Parameter	Value
Product Name	Destrin Rabbit mAb
Remarks/Alias	ADF; ACTDP; HEL32; bA462D18.2; Destrin
Species	Human
GeneID (human)	11034
GeneID	11034
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60981
Protein Weight	19kDa
Shipping	Ice bag