


Dhfr Rabbit Mab For Wb, Elisa - P00374

Item Number: CM0011415



Introduction

DHFR Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat DHFR (21kDa). Suitable for folate metabolism, cancer research, and drug resistance studies.

[Learn More](#)

Parameter	Value
Product Name	DHFR Rabbit mAb
Remarks/Alias	DYR; DHFR1; DHFRP1; DHFR
Species	Human
Gene ID (Human)	1719
Gene ID	1719
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 88-187 of human DHFR (NP_000782.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00374
Protein Weight	21kDa
Shipping	Ice bag