


Dhrs4 Rabbit Pab - Q9Btz2

Item Number: CM0018663



Introduction

High-quality DHRS4 Rabbit Polyclonal Antibody for Western blot and ELISA. Recognizes human, mouse, and rat DHRS4 (~30 kDa), a dehydrogenase/reductase involved in xenobiotic and steroid metabolism. Ideal for IVD assay development.

[Learn More](#)

Parameter	Specification
Product Name	DHRS4 Rabbit pAb
Remarks/Alias	CR; NRDR; PHCR; PSCD; SDR-SRL; SDR25C1; SDR25C2; SCAD-SRL; DHRS4
Species	Human
GeneID (Human)	10901
GeneID	10901
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BTZ2
Protein Weight	30kDa
Shipping	Ice bag