

Digoxin Rabbit Polyclonal Antibody - Digoxin

Item Number: CM0019105

Introduction

Rabbit polyclonal antibody targeting Digoxin for Dot Blot (DB) and ELISA applications. Species-independent reactivity ensures broad utility in digoxin monitoring research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Digoxin Rabbit pAb
Immunogen	Chemical compounds
Source	Rabbit
Category	Polyclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag