


Dimethyl-Eef1A-K55 (P68104/Q05639) Rabbit Polyclonal Antibody For Wb And Elisa

Item Number: CM0000231



Introduction

Rabbit polyclonal antibody recognizing dimethylated lysine 55 of human eEF1A (EEF1A1/EEF1A2, Swiss-Prot P68104/Q05639). Validated for WB and ELISA, ideal for translation elongation factor research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-eEF1A-K55 Rabbit pAb
GeneID (human)	1915/1917
Immunogen	Synthetic dimethylated peptide around K55 of human eEF1A (NP_001393.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P68104/Q05639
Shipping	Ice bag