

Dimethyl-Histone H3-K14 Rabbit Pab (Q16695/P68431) For Wb, If/Icc, Elisa

Item Number: CM0000279



Introduction

Rabbit polyclonal antibody specific for dimethylated lysine 14 of human histone H3 (Q16695/P68431). Suitable for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat, and other species (predicted).

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K14 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K14
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic dimethylated peptide around K14 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag