

Dimethyl-Histone H3-K27 Rabbit Polyclonal Antibody (Q16695/P68431) For Epigenetic Research & Chip

Item Number: CM0000339



Introduction

Rabbit polyclonal antibody specifically recognizing di-methylated lysine 27 of histone H3 (H3K27me2). Suitable for WB, IHC, IF, ELISA, ChIP, and ChIP-seq across human, mouse, and rat. Ideal for epigenetic studies and chromatin dynamics research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K27 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K27
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag