


Dimethyl-Histone H3-K36 (Q16695 / P68431) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Chip-Seq

Item Number: CM0000324



Introduction

Specific rabbit polyclonal antibody recognizing dimethylated Lys36 of human histone H3 (Q16695/P68431). Validated for WB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq applications. Broad cross-reactivity with human, mouse, rat, and other species. Suitable for epigenetic and chromatin biology research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K36 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K36
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic dimethylated peptide around K36 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag