


Dimethyl-Histone H3-K36 Rabbit Monoclonal Antibody (Q16695/P68431) For Wb, Ihc, Ip & Elisa

Item Number: CM0000237



Introduction

High-specificity rabbit monoclonal antibody targeting dimethylated Lys36 of human Histone H3 (SWISS: Q16695/P68431). Validated for WB, DB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, rat and widely predicted species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K36 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K36
Species	Human
GeneID (Human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic dimethylated peptide around K36 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag