


Dimethyl-Histone H3-K4 (Q16695/P68431) Rabbit Monoclonal Antibody For Chip, Wb, And Epigenetic Research

Item Number: CM0000378



Introduction

DiMethyl-Histone H3-K4 rabbit monoclonal antibody targeting human Histone H3 (SWISS Q16695/P68431). Validated for WB, DB, IHC-P, IF/ICC, IP, ELISA, ChIP, ChIP-seq, and CUT&Tag. Cross-reacts with Human, Mouse, Rat, and other species. Suitable for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K4 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K4
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic dimethylated peptide around K4 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, IP, ELISA, ChIP, ChIP-seq, CUT&Tag
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag