


Dimethyl-Histone H3-K4 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Chip, Chip-Seq, Ihc-P

Item Number: CM0000401



Introduction

Rabbit polyclonal antibody targeting dimethylated histone H3 lysine 4 (H3K4me2). Validated for WB, IHC-P, IF/ICC, ELISA, ChIP, and ChIP-seq. Cross-reacts with human, mouse, rat, and other species. Ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K4 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K4
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1))
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag