

Dimethyl-Histone H3-K79 Rabbit Polyclonal Antibody (Q16695/P68431) - Validated For Wb, Ihc-P, Chip-Seq

Item Number: CM0000393



Introduction

Rabbit polyclonal antibody specific for dimethylated Histone H3 at Lysine 79 (H3K79me2). Recognizes human, mouse, rat H3. Suitable for WB, IHC-P, IF, IP, ELISA, ChIP, ChIP-seq. Ideal for epigenetic research and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K79 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K79
Species	Human
GeneID (human)	H3C1, H3C2, H3C3, H3C4, H3C6, H3C7, H3C8, H3C10, H3C11, H3C12, H3-4
GeneID	8290, 8350
Immunogen	A synthetic dimethylated peptide around K79 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag